

24/06/21 - Rep 1 IP

→ Old mouse (w/ lots of fluid and femur bones) + young mouse

→ Bones were crushed thoroughly 6x 5 mL w/ Erylis (RT)

→ 10min RT, filled w/ PBS-only!

→ Spin down 5min, 500xg, 4°C

→ Resuspended in 800µl PBS-only and transferred to low-binding

→ spin down AND resuspended in TRIS lysis buffer: 1 mL

• prepared w/ MS grade buffers (All but 10% NP-40)  
and filtered

• thoroughly vortex until All was sol. and a huge DNA  
pellet was seen

→ spin down 17Kxg, 15min, 4°C

→ sup was collected

→ sup was use w/ 45µl PBS w/ 660nm for  
protein concentration:

y → 2.9 mg/µl

c → 3.9 mg/µl

⇒ IP: will respect the 1mg input: 5µl bead slurry ratio and scale IP reaction accordingly.

y (2.9 mg)  
per IP  
↓  
1.45 mg

α N101 (500µl) + 7.25µl beads  
α IgG (500µl) + 7.25µl beads

o (3.9 mg)  
per IP  
↓  
1.95

α N101 (500µl) + 7.25µl beads  
α IgG (500µl) + 7.25µl beads

- beads were added to new eppies and wash 3x w/ lysis buffer
- IP was incubated 30', rotator, 4°C

⇒ samples were washed 3x 1 ml lysis buffer:

- second wash changed eppies

⇒ After last wash extra spin to remove left over

⇒ Add 25µl 2% SDS 40mM Tris as lysis buffer

⇒ 5 min, 70°C, 1000 rpm

⇒ collect volume in new eppi and store -20°C

06/01/2021 - Rep #2

• 1 young + 1 old (both from Janvier) looking good

→ bones were cleaned and crushed 6x 5ml w/ lysis buffer

• Ellen did the young one

→ lysis 10min, RT, spindown

→ Resuspend in 800µl PBS only

→ transferred to low-binding eppie, spun down

→ y/o samples were resuspended in 1ml lysis buffer  
+  
PI

→ spindown 15 min, 4°C, 17K xg

→ 5µl for 660nm quant:

0.453 mg

Young  $\Rightarrow 453.2529 \mu\text{g/ml} \times 10 = 4.53 \text{ mg/ml}$

Old  $\Rightarrow 460.66134 \mu\text{g/ml} \times 10 = 4.60 \text{ mg/ml}$

0.460 mg

IF 1mg - 5µl

4.6mg - (23µl)

use 23µl bead slurry of αNuc1  
α IgG

- was 23µl beads 3x 1ml RAI's buffer
- there was not enough xNvol so for old was taken from 7 eppie but still same batch.
- 500µl of lysate was added to beads and resuspended
- rotator @ 20, 30min, 4°C
- sort spin down to bring bead down
- 3x 1ml washes w/ eppie change in #2
- extra spin to remove all buffer
- add 25µl 2% SDS 40min RAI's elution buffer
- resuspend beads in elution buffer
- incubate 5min, 70°C, 1000 rpm
- collect elution in new eppie store -20°C

extraction buffer was freshly prepared w/ normal recipe:

- ask leiba about her notes for #1

13/07/2021 - Rep #3

1 young + 1 old (~20mo) from Janvier looking good.

- bones were collected in PBS/BSA
- crushed 6x 5min RT lysis  $\Rightarrow$  Ellen crushed ping one
- 10min RT
- spindown 5min 500xg 4°C
- transferred to low binding eppie
- spin 17K xg, 15min, 4°C
- sup was collected in new eppie
- 5 $\mu$ l of each sample was used for 660nm quant

$$0.46222 \text{ mg} \quad \times 10$$

$$\text{young} \Rightarrow 462.22396 \mu\text{g/ml} \Rightarrow 4.62 \text{ mg/ml} \quad \frac{23.1 \mu\text{l beads}}{\text{add } 23 \mu\text{l}}$$

$$\text{old} \Rightarrow 418.86679 \mu\text{g/ml} \quad \hat{=} 4.18 \text{ mg/ml} \quad \frac{20.9 \mu\text{l beads}}{\text{add } 21 \mu\text{l}}$$

- ⇒ beads were washed 3x 1ml of Tris buffer
- ⇒ 500µl of cell lysate was added to 1P eppie and vortexed
- ⇒ incubation @ 20, 30min, 4°C
- ⇒ beads were washed again 3x 1ml buffer
  - second wash there was eppie change
- ⇒ ~~25~~ extra spin down after first wash to remove excess buf.
- ⇒ 25µl 2% SDS 40mM Tris pH 8.0
- ⇒ 70°C, 5min, 1000 rpm
- ⇒ collect elution into new eppie